

CC152 rabbit pAb

Cat No.:ES19726

For research use only

Overview

Product Name	CC152 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CC152 AA range: 115-165
Specificity	This antibody detects endogenous levels of human CC152
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃. Avoid repeated freeze-thaw cycles.
Protein Name	Coiled-coil domain-containing protein 152
Gene Name	CCDC152 Chr5_400
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	100129792
Human Swiss-Prot Number	Q4G0S7
Alternative Names	
Background	

