



# TRI17 rabbit pAb

Cat No.:ES19734

For research use only

## Overview

Product Name	TRI17 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TRI17 AA range: 95-145
Specificity	This antibody detects endogenous levels of human TRI17
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase TRIM17 (EC 6.3.2.-) (RING finger protein 16) (Testis RING finger protein) (Tripartite motif-containing protein 17)
Gene Name	TRIM17 RBCC RNF16 TERF
Cellular localization	Cytoplasm . Lysosome .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	51127
Human Swiss-Prot Number	Q9Y577
Alternative Names	
Background	

