



# VPS28 rabbit pAb

Cat No.:ES19735

For research use only

## Overview

<b>Product Name</b>	VPS28 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB
<b>Species Cross-Reactivity</b>	Human;Mouse;Rat
<b>Recommended dilutions</b>	WB 1:1000-2000
<b>Immunogen</b>	Synthesized peptide derived from human VPS28 AA range: 19-69
<b>Specificity</b>	This antibody detects endogenous levels of human VPS28
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Vacuolar protein sorting-associated protein 28 homolog (H-Vps28) (ESCRT-I complex subunit VPS28)
<b>Gene Name</b>	VPS28
<b>Cellular localization</b>	Cell membrane . Late endosome membrane ; Peripheral membrane protein .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	
<b>Human Gene ID</b>	51160
<b>Human Swiss-Prot Number</b>	Q9UK41
<b>Alternative Names</b>	
<b>Background</b>	

