

PGPIL rabbit pAb

Cat No.:ES19738

For research use only

Overview

Product Name	PGPIL rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PGPIL AA range: 115-165
Specificity	This antibody detects endogenous levels of human PGPIL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Pyroglutamyl-peptidase 1-like protein (EC 3.4.19.-)
Gene Name	PGPEP1L
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	145814
Human Swiss-Prot Number	A6NFU8
Alternative Names	
Background	

