



EIF2D rabbit pAb

Cat No.:ES19740

For research use only

Overview

Product Name	EIF2D rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human EIF2D AA range: 247-297
Specificity	This antibody detects endogenous levels of human EIF2D
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Eukaryotic translation initiation factor 2D (eIF2d) (Hepatocellular carcinoma-associated antigen 56) (Ligatin)
Gene Name	EIF2D HCA56 LGTN
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	1939
Human Swiss-Prot Number	P41214
Alternative Names	
Background	

