

TPC2A rabbit pAb

Cat No.:ES19741

For research use only

Overview

Product Name	TPC2A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TPC2A AA range: 45-95
Specificity	This antibody detects endogenous levels of human TPC2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Trafficking protein particle complex subunit 2 (Sedlin)
Gene Name	TRAPPC2 SEDL
Cellular localization	Cytoplasm, perinuclear region . Endoplasmic reticulum-Golgi intermediate compartment . Nucleus . Cytoplasm . Localized in perinuclear granular structures. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	6399
Human Swiss-Prot Number	P0DI81
Alternative Names	
Background	

