



STAR6 rabbit pAb

Cat No.:ES19742

For research use only

Overview

Product Name	STAR6 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human STAR6 AA range: 71-121
Specificity	This antibody detects endogenous levels of human STAR6
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	StAR-related lipid transfer protein 6 (START domain-containing protein 6) (StARD6)
Gene Name	STARD6
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	147323
Human Swiss-Prot Number	P59095
Alternative Names	
Background	

