

# MPZL3 rabbit pAb

Cat No.:ES19743

For research use only

## Overview

Product Name	MPZL3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MPZL3 AA range: 59-109
Specificity	This antibody detects endogenous levels of human MPZL3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Myelin protein zero-like protein 3
Gene Name	MPZL3 UNQ2966/PRO7425
Cellular localization	Membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	196264
Human Swiss-Prot Number	Q6UWV2
Alternative Names	
Background	

