



AN36A rabbit pAb

Cat No.:ES19744

For research use only

Overview

Product Name	AN36A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human AN36A AA range: 818-868
Specificity	This antibody detects endogenous levels of human AN36A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Ankyrin repeat domain-containing protein 36A
Gene Name	ANKRD36 ANKRD36A KIAA1641 UNQ2430/PRO499
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	375248
Human Swiss-Prot Number	A6QL64
Alternative Names	
Background	

