

# KLH31 rabbit pAb

Cat No.:ES19745

For research use only

## Overview

Product Name	KLH31 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KLH31 AA range: 342-392
Specificity	This antibody detects endogenous levels of human KLH31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Kelch-like protein 31 (BTB and kelch domain-containing protein 6) (Kelch repeat and BTB domain-containing protein 1) (Kelch-like protein KLHL)
Gene Name	KLHL31 BKLHD6 KBTBD1 KLHL
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	401265
Human Swiss-Prot Number	Q9H511
Alternative Names	
Background	

