



GRP4 rabbit pAb

Cat No.:ES19749

For research use only

Overview

Product Name	GRP4 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GRP4 AA range: 487-537
Specificity	This antibody detects endogenous levels of human GRP4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	RAS guanyl-releasing protein 4
Gene Name	RASGRP4
Cellular localization	Cytoplasm. Cell membrane. Recruited to membranes upon activation by DAG.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	115727
Human Swiss-Prot Number	Q8TDF6
Alternative Names	
Background	

