

ITLN2 rabbit pAb

Cat No.:ES19750

For research use only

Overview

Product Name	ITLN2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ITLN2 AA range: 114-164
Specificity	This antibody detects endogenous levels of human ITLN2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Intelectin-2 (Endothelial lectin HL-2)
Gene Name	ITLN2 UNQ2789/PRO7179
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	142683
Human Swiss-Prot Number	Q8WWU7
Alternative Names	
Background	

