



# TPC3L rabbit pAb

Cat No.:ES19752

For research use only

## Overview

Product Name	TPC3L rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TPC3L AA range: 95-145
Specificity	This antibody detects endogenous levels of human TPC3L
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Trafficking protein particle complex subunit 3-like protein (TRAPPC3-like protein) (BET3-like protein)
Gene Name	BET3L
Cellular localization	Golgi apparatus, cis-Golgi network . Endoplasmic reticulum .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	100128327
Human Swiss-Prot Number	Q5T215
Alternative Names	
Background	

