



PFD3 rabbit pAb

Cat No.:ES19754

For research use only

Overview

Product Name	PFD3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PFD3 AA range: 71-121
Specificity	This antibody detects endogenous levels of human PFD3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Prefoldin subunit 3 (HIBBJ46) (Von Hippel-Lindau-binding protein 1) (VBP-1) (VHL-binding protein 1)
Gene Name	VBP1 PFDN3
Cellular localization	Cytoplasm. Nucleus. In complex with VHL can translocate to the nucleus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7411
Human Swiss-Prot Number	P61758
Alternative Names	
Background	

