

# PCY1B rabbit pAb

Cat No.:ES19757

For research use only

## Overview

Product Name	PCY1B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PCY1B AA range: 89-139
Specificity	This antibody detects endogenous levels of human PCY1B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Choline-phosphate cytidyltransferase B (EC 2.7.7.15) (CCT-beta) (CTP:phosphocholine cytidyltransferase B) (CCT B) (CT B) (Phosphorylcholine transferase B)
Gene Name	PCYT1B CCTB
Cellular localization	[Isoform 1]: Cytoplasm . Endoplasmic reticulum . ; [Isoform 2]: Endoplasmic reticulum .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9468
Human Swiss-Prot Number	Q9Y5K3
Alternative Names	
Background	

