

# ZCHC3 rabbit pAb

Cat No.:ES19761

For research use only

## Overview

Product Name	ZCHC3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZCHC3 AA range: 175-225
Specificity	This antibody detects endogenous levels of human ZCHC3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger CCHC domain-containing protein 3
Gene Name	ZCCHC3 C20orf99
Cellular localization	Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85364
Human Swiss-Prot Number	Q9NUD5
Alternative Names	
Background	

