



DJC16 rabbit pAb

Cat No.:ES19769

For research use only

Overview

Product Name	DJC16 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DJC16 AA range: 171-221
Specificity	This antibody detects endogenous levels of human DJC16
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	DnaJ homolog subfamily C member 16
Gene Name	DNAJC16 KIAA0962
Cellular localization	Membrane ; Single-pass type IV membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23341
Human Swiss-Prot Number	Q9Y2G8
Alternative Names	
Background	

