



MYH7 rabbit pAb

Cat No.:ES19771

For research use only

Overview

Product Name	MYH7 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MYH7 AA range: 1114-1164
Specificity	This antibody detects endogenous levels of human MYH7
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Myosin-7 (Myosin heavy chain 7) (Myosin heavy chain slow isoform) (MyHC-slow) (Myosin heavy chain, cardiac muscle beta isoform) (MyHC-beta)
Gene Name	MYH7 MYHCB
Cellular localization	Cytoplasm, myofibril . Cytoplasm, myofibril, sarcomere . Thick filaments of the myofibrils. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	4625
Human Swiss-Prot Number	P12883
Alternative Names	
Background	

