



# T53G5 rabbit pAb

Cat No.:ES19774

For research use only

## Overview

Product Name	T53G5 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human T53G5 AA range: 85-135
Specificity	This antibody detects endogenous levels of human T53G5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	TP53-target gene 5 protein (TP53-inducible gene 5 protein)
Gene Name	TP53TG5 C20orf10
Cellular localization	Cytoplasm . Nucleus . Cell cycle dependent intracellular localization.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	27296
Human Swiss-Prot Number	Q9Y2B4
Alternative Names	
Background	

