



KRA43 rabbit pAb

Cat No.:ES19777

For research use only

Overview

Product Name	KRA43 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KRA43 AA range: 76-126
Specificity	This antibody detects endogenous levels of human KRA43
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Keratin-associated protein 4-3 (Keratin-associated protein 4.3) (Ultrahigh sulfur keratin-associated protein 4.3)
Gene Name	KRTAP4-3 KAP4.3 KRTAP4.3
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	85290
Human Swiss-Prot Number	Q9BYR4
Alternative Names	
Background	

