

# FOXK2 rabbit pAb

Cat No.:ES19781

For research use only

## Overview

Product Name	FOXK2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FOXK2 AA range: 331-381
Specificity	This antibody detects endogenous levels of human FOXK2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Forkhead box protein K2 (Cellular transcription factor ILF-1) (FOXK1) (Interleukin enhancer-binding factor 1)
Gene Name	FOXK2 ILF ILF1
Cellular localization	Nucleus . Cytoplasm .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3607
Human Swiss-Prot Number	Q01167
Alternative Names	
Background	

