



MAFIP rabbit pAb

Cat No.:ES19789

For research use only

Overview

Product Name	MAFIP rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MAFIP AA range: 48-98
Specificity	This antibody detects endogenous levels of human MAFIP
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	MAFF-interacting protein (Tekin-4 like protein PP5644)
Gene Name	MAFIP MIP PP5644
Cellular localization	Cytoplasm . Nucleus, nucleolus . Interaction with MAFF induces its translocation from the cytoplasm to the nucleolus.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q8WZ33
Alternative Names	
Background	

