

K1161 rabbit pAb

Cat No.:ES19793

For research use only

Overview

Product Name	K1161 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human K1161 AA range: 199-249
Specificity	This antibody detects endogenous levels of human K1161
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Uncharacterized family 31 glucosidase KIAA1161 (EC 3.2.1.-)
Gene Name	KIAA1161
Cellular localization	Nucleus membrane ; Single-pass type II membrane protein . Endoplasmic reticulum membrane ; Single-pass type II membrane protein . Only a minor fraction is present in the peripheral endoplasmic reticulum. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	57462
Human Swiss-Prot Number	Q6NSJ0
Alternative Names	
Background	

