



FBX31 rabbit pAb

Cat No.:ES19794

For research use only

Overview

Product Name	FBX31 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human FBX31 AA range: 17-67
Specificity	This antibody detects endogenous levels of human FBX31
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	F-box only protein 31
Gene Name	FBXO31 FBX14 FBX31 PP2386
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	79791
Human Swiss-Prot Number	Q5XUX0
Alternative Names	
Background	

