

LAP4A rabbit pAb

Cat No.:ES19803

For research use only

Overview

Product Name	LAP4A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LAP4A AA range: 131-181
Specificity	This antibody detects endogenous levels of human LAP4A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lysosomal-associated transmembrane protein 4A (Golgi 4-transmembrane-spanning transporter MTP)
Gene Name	LAPTM4A KIAA0108 LAPTM4 MBNT MTRP UNQ1846/PRO3574
Cellular localization	Endomembrane system ; Multi-pass membrane protein . May reside in an intracellular membrane-bound compartment. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9741
Human Swiss-Prot Number	Q15012
Alternative Names	
Background	

