



CK037 rabbit pAb

Cat No.:ES19808

For research use only

Overview

Product Name	CK037 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human CK037 AA range: 55-105
Specificity	This antibody detects endogenous levels of human CK037
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Putative uncharacterized protein encoded by LINC00167
Gene Name	LINC00167 C11orf37 NCRNA00167
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q96N53
Alternative Names	
Background	

