



PNMA2 rabbit pAb

Cat No.:ES19813

For research use only

Overview

Product Name	PNMA2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human PNMA2 AA range: 153-203
Specificity	This antibody detects endogenous levels of human PNMA2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Paraneoplastic antigen Ma2 (40 kDa neuronal protein) (Onconeuroal antigen Ma2) (Paraneoplastic neuronal antigen MM2)
Gene Name	PNMA2 KIAA0883 MA2
Cellular localization	Nucleus, nucleolus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10687
Human Swiss-Prot Number	Q9UL42
Alternative Names	
Background	

