



GBRL1 rabbit pAb

Cat No.:ES19814

For research use only

Overview

Product Name	GBRL1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GBRL1 AA range: 32-82
Specificity	This antibody detects endogenous levels of human GBRL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Gamma-aminobutyric acid receptor-associated protein-like 1 (Early estrogen-regulated protein) (GABA(A) receptor-associated protein-like 1) (Glandular epithelial cell protein 1) (GEC-1)
Gene Name	GABARAPL1 GEC1
Cellular localization	Cytoplasmic vesicle, autophagosome . Cytoplasmic vesicle membrane ; Lipid-anchor . Cytoplasm, cytoskeleton . Endoplasmic reticulum . Golgi apparatus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	23710
Human Swiss-Prot Number	Q9H0R8
Alternative Names	
Background	

