

SNX20 rabbit pAb

Cat No.:ES19816

For research use only

Overview

Product Name	SNX20 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SNX20 AA range: 90-140
Specificity	This antibody detects endogenous levels of human SNX20
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sorting nexin-20 (Selectin ligand-interactor cytoplasmic 1)
Gene Name	SNX20 SLIC1
Cellular localization	Early endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Cell membrane . Cytoplasm . Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	124460
Human Swiss-Prot Number	Q7Z614
Alternative Names	
Background	

