

HDX rabbit pAb

Cat No.:ES19817

For research use only

Overview

Product Name	HDX rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human HDX AA range: 321-371
Specificity	This antibody detects endogenous levels of human HDX
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Highly divergent homeobox
Gene Name	HDX CXorf43
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	139324
Human Swiss-Prot Number	Q7Z353
Alternative Names	
Background	

