

TM86A rabbit pAb

Cat No.:ES19818

For research use only

Overview

Product Name	TM86A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TM86A AA range: 154-204
Specificity	This antibody detects endogenous levels of human TM86A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Lysoplasmalogenase-like protein TMEM86A (Transmembrane protein 86A)
Gene Name	TMEM86A
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	144110
Human Swiss-Prot Number	Q8N2M4
Alternative Names	
Background	

