



TMEM9 rabbit pAb

Cat No.:ES19819

For research use only

Overview

Product Name	TMEM9 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TMEM9 AA range: 29-79
Specificity	This antibody detects endogenous levels of human TMEM9
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Transmembrane protein 9 (Dermal papilla-derived protein 4)
Gene Name	TMEM9 DERP4 HSPC186 PSEC0012 UNQ631/PRO1248
Cellular localization	Lysosome membrane ; Single-pass type I membrane protein . Late endosome membrane ; Single-pass type I membrane protein . Endosome, multivesicular body membrane ; Single-pass type I membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	252839
Human Swiss-Prot Number	Q9P0T7
Alternative Names	
Background	

