

TJAP1 rabbit pAb

Cat No.:ES19822

For research use only

Overview

Product Name	TJAP1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human TJAP1 AA range: 447-497
Specificity	This antibody detects endogenous levels of human TJAP1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Tight junction-associated protein 1 (Protein incorporated later into tight junctions) (Tight junction protein 4)
Gene Name	TJAP1 PILT TJP4
Cellular localization	Golgi apparatus, trans-Golgi network . Cell junction, tight junction . Cell membrane ; Peripheral membrane protein . Recruited to tight junctions (TJ) during late stages of maturation of the TJ complexes. Excluded from adherens junctions and desmosomes. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	93643
Human Swiss-Prot Number	Q5JTD0
Alternative Names	
Background	

