



RN215 rabbit pAb

Cat No.:ES19823

For research use only

Overview

Product Name	RN215 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human RN215 AA range: 104-154
Specificity	This antibody detects endogenous levels of human RN215
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	RING finger protein 215
Gene Name	RNF215
Cellular localization	Membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	200312
Human Swiss-Prot Number	Q9Y6U7
Alternative Names	
Background	

