

IVD rabbit pAb

Cat No.:ES19827

For research use only

Overview

Product Name	IVD rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human IVD AA range: 325-375
Specificity	This antibody detects endogenous levels of human IVD
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Isovaleryl-CoA dehydrogenase, mitochondrial (IVD) (EC 1.3.8.4)
Gene Name	IVD
Cellular localization	Mitochondrion matrix .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	3712
Human Swiss-Prot Number	P26440
Alternative Names	
Background	

