

UDB11 rabbit pAb

Cat No.:ES19828

For research use only

Overview

Product Name	UDB11 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human UDB11 AA range: 53-103
Specificity	This antibody detects endogenous levels of human UDB11
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20℃ . Avoid repeated freeze-thaw cycles.
Protein Name	UDP-glucuronosyltransferase 2B11 (UDPGT 2B11) (EC 2.4.1.17)
Gene Name	UGT2B11
Cellular localization	Microsome membrane ; Single-pass membrane protein . Endoplasmic reticulum membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10720
Human Swiss-Prot Number	O75310
Alternative Names	
Background	

