

LRGUK rabbit pAb

Cat No.:ES19831

For research use only

Overview

Product Name	LRGUK rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human LRGUK AA range: 688-738
Specificity	This antibody detects endogenous levels of human LRGUK
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Leucine-rich repeat and guanylate kinase domain-containing protein
Gene Name	LRGUK
Cellular localization	Cytoplasmic vesicle, secretory vesicle, acrosome . Cytoplasm, cytoskeleton . Cytoplasm, cytoskeleton, cilium basal body . Localizes to the acrosome and acroplaxome in round spermatids. Localizes to the manchette during spermiogenesis. Also found in the basal body of elongating spermatids, and in primary cilia of somatic cells. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	136332
Human Swiss-Prot Number	Q96M69
Alternative Names	
Background	

