

SDF2 rabbit pAb

Cat No.: ES19833

For research use only

Overview

Product Name SDF2 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human;Mouse **Recommended dilutions** WB 1:1000-2000

Immunogen Synthesized peptide derived from human SDF2 AA

range: 152-202

Specificity This antibody detects endogenous levels of human

SDF2

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20° C. Avoid repeated freeze-thaw cycles.

Protein Name Stromal cell-derived factor 2 (SDF-2)

Gene Name SDF2 **Cellular localization** Secreted .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

Observed band

Human Gene ID 6388 Human Swiss-Prot Number Q99470

Alternative Names

Background

