

## ZHX1 rabbit pAb

Cat No.: ES19834

For research use only

## Overview

Product Name ZHX1 rabbit pAb

Host species Rabbit
Applications WB

Species Cross-Reactivity Human; Mouse; Rat Recommended dilutions WB 1:1000-2000

Immunogen Synthesized peptide derived from human ZHX1 AA

range: 736-786

**Specificity** This antibody detects endogenous levels of human

ZHX1

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

**Storage** Store at  $-20^{\circ}$ C. Avoid repeated freeze-thaw cycles.

Protein Name Zinc fingers and homeoboxes protein 1

Gene Name ZHX1

**Cellular localization**Nucleus . Colocalized in the nucleus with DNMT3B. **Purification**The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal Concentration 1 mg/ml

**Observed band** 

Human Gene ID 11244 Human Swiss-Prot Number Q9UKY1

**Alternative Names** 

+86-27-59760950

**Background** 

