



ZFPL1 rabbit pAb

Cat No.:ES19843

For research use only

Overview

Product Name	ZFPL1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZFPL1 AA range: 240-290
Specificity	This antibody detects endogenous levels of human ZFPL1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein-like 1 (Zinc finger protein MCG4)
Gene Name	ZFPL1
Cellular localization	Golgi apparatus, cis-Golgi network membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7542
Human Swiss-Prot Number	O95159
Alternative Names	
Background	

