



AP4M1 rabbit pAb

Cat No.:ES19845

For research use only

Overview

Product Name	AP4M1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human AP4M1 AA range: 312-362
Specificity	This antibody detects endogenous levels of human AP4M1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	AP-4 complex subunit mu-1 (AP-4 adapter complex mu subunit) (Adapter-related protein complex 4 mu-1 subunit) (Mu subunit of AP-4) (Mu-adaptin-related protein 2) (mu-ARP2) (Mu4-adaptin) (mu4)
Gene Name	AP4M1 MUARP2
Cellular localization	Golgi apparatus, trans-Golgi network membrane ; Peripheral membrane protein . Early endosome . Found in soma and dendritic shafts of neuronal cells. .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	9179
Human Swiss-Prot Number	O00189
Alternative Names	
Background	

