

NICN1 rabbit pAb

Cat No.:ES19846

For research use only

Overview

Product Name	NICN1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human NICN1 AA range: 145-195
Specificity	This antibody detects endogenous levels of human NICN1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Nicotin-1 (NPCEDRG) (Tubulin polyglutamylase complex subunit 5) (PGs5)
Gene Name	NICN1
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84276
Human Swiss-Prot Number	Q9BSH3
Alternative Names	
Background	

