



DJC30 rabbit pAb

Cat No.:ES19847

For research use only

Overview

Product Name	DJC30 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human DJC30 AA range: 50-100
Specificity	This antibody detects endogenous levels of human DJC30
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	DnaJ homolog subfamily C member 30 (Williams-Beuren syndrome chromosomal region 18 protein)
Gene Name	DNAJC30 WBSCR18
Cellular localization	Mitochondrion inner membrane ; Single-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	84277
Human Swiss-Prot Number	Q96LL9
Alternative Names	
Background	

