



ZNF2 rabbit pAb

Cat No.:ES19849

For research use only

Overview

Product Name	ZNF2 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ZNF2 AA range: 25-75
Specificity	This antibody detects endogenous levels of human ZNF2
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	Zinc finger protein 2 (Zinc finger protein 2.2) (Zinc finger protein 661)
Gene Name	ZNF2 ZNF661
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	7549
Human Swiss-Prot Number	Q9BSG1
Alternative Names	
Background	

