



KIAS3 rabbit pAb

Cat No.:ES19857

For research use only

Overview

Product Name	KIAS3 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human KIAS3 AA range: 62-112
Specificity	This antibody detects endogenous levels of human KIAS3
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Putative uncharacterized protein KIRREL3-AS3 (KIRREL3 antisense RNA 1) (KIRREL3 antisense gene protein 1) (Proline-rich protein 10)
Gene Name	KIRREL3-AS3 NCRNA00288 PRR10
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	
Human Swiss-Prot Number	Q8N7Y1
Alternative Names	
Background	

