

ADRO rabbit pAb

Cat No.:ES19859

For research use only

Overview

Product Name	ADRO rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human ADRO AA range: 364-414
Specificity	This antibody detects endogenous levels of human ADRO
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	NADPH:adrenodoxin oxidoreductase, mitochondrial (AR) (Adrenodoxin reductase) (EC 1.18.1.2) (Ferredoxin--NADP(+) reductase) (Ferredoxin reductase)
Gene Name	FDXR ADXR
Cellular localization	Mitochondrion inner membrane ; Peripheral membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	2232
Human Swiss-Prot Number	P22570
Alternative Names	
Background	

