



SG11B rabbit pAb

Cat No.:ES19860

For research use only

Overview

Product Name	SG11B rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SG11B AA range: 38-88
Specificity	This antibody detects endogenous levels of human SG11B
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Sperm-associated antigen 11B (Human epididymis-specific protein 2) (He2) (Protein EP2) (Sperm antigen HE2)
Gene Name	SPAG11B EP2 HE2
Cellular localization	Secreted .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10407
Human Swiss-Prot Number	Q08648
Alternative Names	
Background	

