

SUGT1 rabbit pAb

Cat No.:ES19862

For research use only

Overview

Product Name	SUGT1 rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human SUGT1 AA range: 96-146
Specificity	This antibody detects endogenous levels of human SUGT1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	Suppressor of G2 allele of SKP1 homolog (Protein 40-6-3) (Sgt1)
Gene Name	SUGT1
Cellular localization	Cytoplasm . Nucleus . Translocates to the nucleus upon heat shock, requiring S100A6.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	10910
Human Swiss-Prot Number	Q9Y2Z0
Alternative Names	
Background	

