

GAG2A rabbit pAb

Cat No.:ES19864

For research use only

Overview

Product Name	GAG2A rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human GAG2A AA range: 48-98
Specificity	This antibody detects endogenous levels of human GAG2A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	G antigen 2A (GAGE-2A)
Gene Name	GAGE2A GAGE2
Cellular localization	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	729447
Human Swiss-Prot Number	Q6NT46
Alternative Names	
Background	

