

MARHB rabbit pAb

Cat No.:ES19872

For research use only

Overview

Product Name	MARHB rabbit pAb
Host species	Rabbit
Applications	WB
Species Cross-Reactivity	Human;Mouse;Rat
Recommended dilutions	WB 1:1000-2000
Immunogen	Synthesized peptide derived from human MARHB AA range: 198-248
Specificity	This antibody detects endogenous levels of human MARHB
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C . Avoid repeated freeze-thaw cycles.
Protein Name	E3 ubiquitin-protein ligase MARCH11 (EC 6.3.2.-) (Membrane-associated RING finger protein 11) (Membrane-associated RING-CH protein XI) (MARCH-XI)
Gene Name	44266
Cellular localization	Cytoplasmic vesicle membrane ; Multi-pass membrane protein .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	
Human Gene ID	441061
Human Swiss-Prot Number	A6NNE9
Alternative Names	
Background	

